


PRODUCT LAUNCH Promo - BUY 1 FREE 1



- ✓ Angled hose connector simplifies vacuum line attachment
- ✓ Easy grip sides receiver bottles for improved handling
- ✓ Individually packaged, sterile
- ✓ Certified nonpyrogenic

Important Performance Characteristics of Corning Filter Membrances



PES=Polyethersulfone, CA=Cellulose Acetate, CN=Cellulose Nitrate, NY=Nylon



1 pack of 0.2µm MiniSart Syringe Filter with ANY 5 cases purchased in a single order

CORNING Vacuum Filter System Ordering Information

Catalog Number	Description	Packing
431096	250ml PES Filter unit, 50mm diameter membrane, 0.22µm pore size	12/ Case
431097	500ml PES Filter unit, 50mm diameter membrane, 0.22µm pore size	12/ Case
431098	1000ml PES Filter unit, 50mm diameter membrane, 0.22µm pore size	12/ Case
431118	500ml PES Bottle Top Filter, 70mm diameter membrane, 45mm neck 0.22µm pore size	12/ Case

* Offer valid while stocks last.