


***SPECIAL OFFER !!!**



Highest Capacity Refrigerated Micro Centrifuge: 3 x 24 x 1.5mL



Rack-in-Rotor

Stackable, drop-in/lift-out rotor racks for easy placement & removal

Lightweight resin racks have a drop-in/lift-out design which allows for easy placement and replacement of tubes. These rotor racks are available in a wide variety of types, and can be stacked in layers to suit your application.

Lightweight, yet stable rotation, and excellent acceleration/ deceleration performance

Each rack has four slits so that the adhesion of the rack to the inner wall of the rotor is facilitated by centrifugal force, providing stable operation of the rotor. Lightweight rotors and racks can also increase acceleration and deceleration performance.

A wide variety of racks available at affordable prices

A large selection of resin racks is provided at affordable prices to accommodate tubes and plates of different volumes and numbers, which supports a diverse array of sample processing applications.

A portable rack for sample tubes

A resin rack can be used for a wide range of purposes, including storage and transportation. It can be used as a tube stand or carrier to hold or transport sample tubes.



TMA-300



AR510-04



AR015-24

Easy-to-see, easy-to-use control panel

During stopping



During spinning



- User friendly center Jog Dial operation
- Highly visible dual color LED display for operating status indication
- Memory function storing up to three sets of operating parameters for more efficient centrifugation
- One touch key operation permits easy acceleration of the rotor to a maximum speed for the centrifugation protocol selected
- 1 - 99 minutes timer with 1 minute increment

Contact our sales representatives
for this special offer ...

* Offer valid while stock last